# REQUEST FOR LEGISLATIVE ACTION

Completed by County Counselor's Office:
Res/Ord No.: 20070
Sponsor(s): Crystal Willia
Date: December 3, 20

Crystal Williams December 3, 2018

SUBJECT	Action Requested Resolution Ordinance	
	Project/Title: Requesting the transfer of \$68,743 within the Medical Examiner's budg equipment.	et for the purchase of digital xray
BUDGET		
INFORMATION	Amount authorized by this legislation this fiscal year:	\$68,743.
To be completed	Amount previously authorized this fiscal year:	\$0
By Requesting	Total amount authorized after this legislative action:	\$68,743.
Department and Finance	Amount budgeted for this item * (including transfers):	\$0
r munce	Source of funding (name of fund) and account code number; FROM: 002-2001 56570 Health Fund, MEO, Maint & Repair Misc	FROM ACCT \$68,743.
	TO: 002-2001 58170 Health Fund, MEO, Other Equipment	TO ACCT \$68,743.
	OTHER FINANCIAL INFORMATION:  No budget impact (no fiscal note required) Term and Supply Contract (funds approved in the annual budget); es Department: Estimated Use: \$  Prior Year Budget (if applicable): \$  Prior Year Actual Amount Spent (if applicable): \$	stimated value and use of contract:
PRIOR LEGISLATION.	Prior ordinances and (date): Prior resolutions and (date): Ord 5174 11/26/18	
CONTACT INFORMATION	RLA drafted by (name, title, & phone): Kandi Brooke, Administrative S Dr. Diane Peterson, Chief Medical Examiner (816) 881-6600	Supervisor for
REQUEST SUMMARY	Requesting the transfer of \$68,743 within the Medical Examiner's budge equipment DR System, GE AMX 4, and Fuji T2 CR Imaging	et for the purchase of digital xray
CLEARANCE	Tax Clearance Completed (Purchasing & Department) Business License Verified (Purchasing & Department) Chapter 6 Compliance - Affirmative Action/Prevailing Wage (Count	ty Auditor's Office)

ATT	ACHMENTS	Quote for DR System (S	\$27,943.75), GE AMX 4 (\$	19,000), and Fuji T2 CR Imaging	(\$21,800)
REV	EW	Department Director:	Drane Ost	teronMD	Date: 11/28/2018
		Finance (Budget Appro If applicable	val):		Date: 11/2018
		Division Manager:	Daine M		Date: 11-28-18
		County Counselor's Off	fice:		Date:
		diture was included in the	Budget Office in Finan		
	Funds for the	nis were encumbered from	n the	Fund in	
₩.	is chargeab	le and there is a cash bala		appropriation to which the expended in the treasury to the credit of ion herein authorized.	
	Funds suffi	cient for this expenditure	will be/were appropriated	by Ordinance #	
	Funds suffi	cient for this appropriation	n are available from the so	urce indicated below.	
	Account N	Jumber:	Account Title:	Amount Not to Exce	ed:
				on County to pay any specific ame	

This legislative action does not impact the County financially and does not require Finance/Budget approval.

# **Fiscal Note:**

Funds sufficient for this transfer are available from the sources indicated below.

	'							
Date:	November 28, 2018					RES#_	2007	0
Depart	ment / Division	Chara	cter/Description		Froi	m	То	
002	Health Fund							
2001	Medical Examiner	56570	Maint. & Repair - M	isc.	\$	68,743	\$	
2001	Medical Examiner	58170	Other Equipment		<del></del>	-	ş <del></del>	68,743
					÷		¥	
					·		¥ <del></del>	
				=:				
							13	
				<del>.</del>		***		
							3	
		). <del>(1</del>			<del>.</del>		9	
		-						
			±1	<del></del>	5 <del></del>		3	
		<del></del>					-	
	11				\$	68,743	\$	68,743

Budget Officer

### OFFICE OF THE JACKSON COUNTY MEDICAL EXAMINER

950 East 21st Street Kansas City, Missouri 64108 (816) 881-6600 (816) 881-6641 fax

To: Barbara Casamento, Purchasing

From: Diane C. Peterson, MD, Chief Medical Examiner

Re: Emergency Purchase of new x-ray equipment

The Medical Examiner Office had a previous ordinance to repair our broken x-ray equipment (Ordinance 5174). However, in the time between initial request and passage, our office was made aware of the limited lifespan of the machine, even after repair. Therefore, we request funds to be moved from 6570 Maintenance & Repair to 8170 Other Equipment for the emergency purchase of new machine and readers. We have been using our KCRMORG x-ray to get by, but it is not ideal as it does not fit our x-ray table.

We would like to complete this purchase this budget year with currently available funds. Transfer of funds within our budget is sufficient. No outside funds will be necessary.

If you have any questions, please let me know.

Thank you,

Diane C. Peterson, MD Chief Medical Examiner



8121 NW Terrace, Kansas City, MO 64153 Phone: 816-741-5558 Fax: 816-587-2423 # 19,000. 00
Visit Us Online At:

Visit Us Online At: www.radsource.net

### **SALES QUOTATION**

QUOTE # RADQ4225

DATE 11/8/2018

VALID THROUGH 12/8/2018



Account Manager Nancy Briney
Phone 816-741-5558
Email nbriney@radsource.net

### **CUSTOMER DETAILS**

### **Jackson County Medical Examine**

Bill To:

JACKSON CO MED EXAM-2062-MO 950 EAST 21ST STREET KANSAS CITY, MO 64108

Dr Diane Peterson / Medical Examin

Ship To:

Phone (816)881-6600

Email KBrooke@jacksongov.org

GPO:

Part Number Unit Price Qty Ext. Price

### Platinum Select GE AMX-4+ Portable System

#### **GE AMX 4 Portable Include:**

- MX-75 275k HU X-ray tube with a .75mm focal spot
- · GE AMX-IV collimator
- · High-quality image
- · Forward and reverse variable speed
- · Completely mobile
- LED display
- · Operates at full power with no power connection
- 24 kVp stations
- · 30 mAs stations
- · 100 mA constant
- kVp 50-125
- · mAs 0.4-320

10,000

Initial Here

#### **Dimensions**

• Size (height x width x length): 76" (1,930 mm) x 25 3/16 (640 mm) x 45 3/8" (1.153 mm)

• Weight: 1,050 lbs (477 kg)

#### **Battery**

Nine 12.9 volt batteries connected in series provide approximately 116 volt at full charge. Battery capacity is approximately 20,000 mAs at 100 kVp when fully charged and under specific conditions. Actual mAs capacity depends on the condition of the batteries, use of collimator lamp, electric brakes, rotor, and distance driven.

#### Movements (tube vertical movement measured at the focal spot)

- Range: At least 52.5" (1,334 mm)
- Lowest position: 26.1" (663 mm) maximum from floor
- Highest position: 78.6" (1,996 mm) minimum from floor
- Horizontal movement measured at focal spot relative to column face: 24" (610 mm) minimum to 40" (1.016 mm) maximum
- Tube column rotation measured from horizontal latch: 360°
- Tube and yoke rotation around horizontal arm measured from tube port down position: Range (360°), detent locations (0, ±90°, ±180°)
- Tube trunion rotation measured from tube port down position: Range (120°), forward (110°), backward (10°), detent (0°, 90°)
- Collimator rotation measured from the front of the collimator with the tube port facing down: Range (180°), right (90°), left (90°), detent (0°, 90°)

#### **Drive Speeds**

Speed 1: Drive speed with horizontal arm secured for support - 264' (6,705 mm) per minute ±25%

Speed 2: Maneuvering speed with horizontal tube arm removed from the transport latch - (30% to 60% of drive speed)

\*Drive speed is measured on a smooth, hard and level surface. Speed will be reduced by inclines, carpeted or soft surfaces.

#### Reproducibility

Coefficient of variation of radiation output is less than 0.045 for successive exposures having constant technique factors.

#### Collimator

Minimum source to skin distance is limited to more than 30 cm by the skin spacers at the beam exit of the collimator. Full 17" x 17" (43 cm x 43 cm) radiographic coverage at 40" (1.02 mm) source to image distance. Minimum inherent filtration of 2.0 mm aluminum equivalent at 100 kVp.

\*Note:If purchased with a Used CR or New DR panel AMX IV will be priced at \$19,000

19000

	System	Purchase Price:	\$20,000.00
Taxes and Freight charges will be included in specified on this quotation.	the final invoice if not	Subtotal Sales Tax	\$20,000.00
Tax Exempt: YES / NO (please circle one) If yes, we will require a copy of your tax exempt certificate.		Shipping Grand Total	\$20,000.00
Please choose / check by one of the following	g:	Deposit Required	\$10,000.00
Funding Source:  Cash Credit Card Bank Loan Lease	Credit Card Name:		order is required.

#### **Payment Terms:**

The terms stated in this agreement are an important part of this business/purchase agreement. The buyer certifies and agrees to the stated payment terms and agrees that they or their purchasing organization will carry them out as stated in this agreement.

50% due with signed buyer's agreement, 40% due 10-Business Days prior to the scheduled installation of this equipment, remaining 10% due upon completion of installation.

Additional convenience fee of 3% will be added to all invoices paid by approved credit card.

Pricing contained within this proposal excludes shipping and applicable taxes. Shipping charges are pre-paid by RadSource and shall be added to the buyer's final invoice. Any applicable taxes are the responsibility of the buyer. If the buyer is tax-exempt, RadSource Imaging will require a copy of their tax exempt certificate.

System(s) installation will be scheduled on a mutually agreeable date. Installation will not begin until the (first and second) progress payments have been submitted to RadSource Imaging Technologies, Inc.

### Warranty:

GE AMX IV PARTS WARRANTY:

Twelve (36) months full warranty, for all parts, including X-ray tube, and batteries.

All warranties are subject to manufacturer's terms and conditions. Service hours of warranty coverage are Monday through Friday from 8:00 AM to 5:00 PM excluding holidays.

Extended warranty Information available under separate cover.

#### **Warranty Exclusions:**

DR Drop Protection excluded if applicable. Warranty does not cover damage caused by any act of God, war, terrorism, normal wear and tear, dropping or impact, neglect, power spikes or surges that might potentially damage the equipment's internal components, loss of power, water, smoke, fire, fire suppressant, high winds, lightning strikes, viruses, tampering, operator abuse or use other than what the equipment is intended to be used for. Warranty is non-transferrable. RadSource Imaging Technologies, Inc. is not responsible for lost revenue due to down time.

DR panel drop protection not included in system warranty.

#### Notes:

Upon acceptance by both parties, as indicated by their signatures below, this Quotation shall become a binding Sales Contract between BUYER and SELLER.

Installation will be performed during RadSource Imaging Technologies normal operational hours, Monday through Friday, 8:00AM to 5:00PM. Requests for overtime installation srevices will need to be approved by RadSource Imaging Technologies Service Manager. RadSource Imaging Technologies prevailing overtime rates will apply, times the number of men required to perform the required installation.

#### **Buyer Site Preparation and Installation:**

Buyer is responsible for all site and/or room preparation including but not limited to Electrical, Plumbing, Overhead Support Work, Power Access, and providing the necessary Networking Access. Professional room drawings can be provided by RadSource Imaging Technologies Inc. upon customer's request with purchase of equipment from RadSource. Drawings will be started once down payment is received. The area where the X-ray equipment is to be installed will need to be completely finished including paint, baseboards, ceiling etc. before the equipment installation can take place.


Quoted Prices and specifications are subject to change without notice prior to signed buyers agreement. Service Agreement Options are available to protect your purchase after the initial warranty has expired.

During or following installation, a RadSource employee may take photographs of the system being installed at your facility along with an image of the outside of your facility to include on their website's installation gallery

(www.radsource.net/installation-gallery/). By signing this agreement, you agree to grant to RadSource, permission to use, reuse, publish, and republish any photographs that they have taken, or may take of your facility, or the equipment installed at your facility. (RadSource will never take photographs of people, or any other hospital, patient or other confidential information.) The photographs can be used in any and all media now or hereafter, including the internet, for the purpose of promoting our installation process at your facility; and to use your company name in connection with any images RadSource so chooses to use. You agree to release and discharge RadSource from any and all claims and demands that may arise out of or in connection with the use of the photographs.

CUSTOMER ACCEPTANCE	RADSOURCE ACCEPTANCE
NAME	NAME
SIGN	SIGN
DATE	DATE

Thank You For Your Business!





8121 NW Terrace, Kansas City, MO 64153 Phone: 816-741-5558 Fax: 816-587-2423 Visit Us Online At: www.radsource.net

### **SALES QUOTATION**

QUOTE # RADQ4226

DATE 11/8/2018

VALID THROUGH 12/8/2018



Account Manager Nancy Briney
Phone 816-741-5558
Email nbriney@radsource.net

### **CUSTOMER DETAILS**

## **Jackson County Medical Examine**

Bill To:

JACKSON CO MED EXAM-2062-MO 950 EAST 21ST STREET KANSAS CITY, MO 64108

Dr Diane Peterson / Medical Examin

Ship To:

Phone (816)881-6600

Email KBrooke@jacksongov.org

GPO:

Part Number Unit Price Qty Ext. Price

New Vieworks VIVIX-S DR 14x17V Wireless Gadox DR System

HV6056

VIVIX .S View DR 14x17 Wireless Gadox Panel . Scintillator: Amorphous Silicon Gadolinium Oxysulfide Gd.O.S:Tb (Gadox) \$27,200.00

1

\$27,200.00

Scintillator Amorphous Silicon: Gadolinium Oxysulfide Gd.O.S:Tb (Gadox)

ALTO DR™ 14x17 Wireless Detectors Include: 14x17 Wireless Flat Panel Detector, System Control Unit (SCU), Battery Charger, Battery Pack, Image Processing Viewer Software, Cables, and Acquisition Workstation.

Detection Area: 14x17 inch (35cm x 43cm)
Dimensions: 18.1"W x 15.1"L x .59"D (46 x 38.4 x 1.5 cm)

- Pixels: 2560 x 3072
- Pixel Pitch: 140 im x 140 im
- Scintillator: Cesium (CsI:TI) or Gadox (Gd2O2S:Tb)
- · Grayscale: 16,384
- A/D Conversion: 14 bit
- X-Ray Voltage Range: 40-150 kVp
- Weight: 7.27 lbs./3.3 kg. (CsI) or 6.83 lbs./3.1 kg. (Gadox)
- Power: 7.4V 4000mAh Lithium Ion Polymer Battery

DR ACQUISION WORKSTATION

### **Acquisition Workstation**

Acquisition Hardware:

- Intel Core i3 Processor 6100 CPU (3.7GHz)
- 4 GB. Windows 7 Professional 64-bit
- 1 TB Local Mini-PACs in a RAID 1 (2 drives total)
- Professional Monitor: 22 in. Touch Screen Resolution: 1920 X 1080, full HD Ratio: 1000:1 (typical)

VxVue Imaging Acquisition Software:

- Viewers 5 (LAN) Local Area Network Clinical Viewing Licenses
- AED Technology Automatic Exposure Detection: No Generator Interfacing is required, as the AED Feature instantly detects x-ray upon exposure and then completes the acquisition process for review and transfer of studies to PACS
- DICOM Store/Transfer to PACS, Auto-Routing and Scheduled periodic back-up
- Image Preview in 3 seconds
- Customizable Anatomical exam tree and image processing parameters.
- Image Acquisition + QC (brightness/contrast, window/level, flip, rotate, zoom, annotations, measurements, pan, mag, crop, invert, line measurement, angle measurement, text, arrow, L&R Markers, clockwise rotation, counter-clockwise rotation.)
- Patient Registration Manually or DICOM Modality Worklist. DICOM Modality Worklist is a standard feature and is available for use provided that the EMR/EHR/RIS accommodates Modality Worklist.
- Print: DICOM Printer or Desk Top Paper Printer.
- Stitching: Manual and Auto.
- Archive: CD/DVD Archive to removable drive.
- Export Images as, JPEG, TIFF, BMP, and DICOM with or without Viewer to CD or other media.
- Patient Size Selection: Pediatric, Small, Medium, Large Adult.
- Black Surround Fill

APPS 2.0 2 Days Applications Training Included in purchase price.

**System Purchase Price** 

\$27,200.00

Part Number		Unit Price	Qty	Ext. Price	
	Optional Items (Please add to system pu	irchase price	)		
WDDR14	14X17 WEIGHT DISTRIBUTION CAP (Optional)	\$400.00	) 1	Initials	
PH1417	14X17 PROTECT-A-GRID (PLASTIC) ENCASEMENT WITH 2 HANDLES AND DOUBLE RELEASE LOCKS (Optional)	\$743.75	) 1 .	PP Initials	
K14102040SD	Must be ordered with the below "Installed Grid" Cat. #K1410 Please add both prices for the complete Protect-A-Grid Protect Appropriate Grid set. Specified Grid is built into a durable pla protective casing.  14X17 10:1/215/40-72 GRID (SHORT DIMENSION) (Optional)	ctor and	7	Initials	
TTDR	TABLETOP DR PANEL HOLDER (Optional) \$600.00 1 Initial				
	(used for performing lateral tabletop radiographs. Holds the  Protect-A-Grid and DR Panel)				
	If DR is purchased with AMX IV + from RadSource the DR price will be lowered to \$26,200				
Taxes and Fre	ight charges will be included in the final invoice if not	Subtotal	- X - 3	\$27,200.00	
specified on th		Sales Tax			
Tax Exempt:	YES / NO (please circle one)	Shipping			
	require a copy of your tax exempt certificate.	Grand Total		\$27,200.00	
Please choose	/ check by one of the following:	Deposit Rec	uired	\$0.00	
Funding Source	e: Cash		(2	7943.7	

\*If order is being leased, a Leasing Company's Purchase Order is required.

### **Payment Terms:**

The terms stated in this agreement are an important part of this business/purchase agreement. The buyer certifies and agrees to the stated payment terms and agrees that they or their purchasing organization will carry them out as stated in this agreement.

50% due with signed buyer's agreement, 40% due 10-Business Days prior to the scheduled installation of this equipment, remaining 10% due upon completion of installation.

Lessor:

Additional convenience fee of 3% will be added to all invoices paid by approved credit card.

Lease

Pricing contained within this proposal excludes shipping and applicable taxes. Shipping charges are pre-paid by RadSource and shall be added to the buyer's final invoice. Any applicable taxes are the responsibility of the buyer. If the buyer is tax-exempt, RadSource Imaging will require a copy of their tax exempt certificate.

System(s) installation will be scheduled on a mutually agreeable date. Installation will not begin until the (first and second) progress payments have been submitted to RadSource Imaging Technologies, Inc.

### Warranty:

### VIEWORKS ALTO MANUFACTURERS WARRANTY:

Sixty (60) Months on manufacturer defects on the DR Detector. 36-Months on manufacturer defects on the standard PC acquisition station, and software.

RADQ4226

Initial Here

#### RADSOURCE LABOR WARRANTY:

Twelve (12) months service from date of completion of installation.

All warranties are subject to manufacturer's terms and conditions. Service hours of warranty coverage are Monday through Friday from 8:00 AM to 5:00 PM excluding holidays.

Extended warranty Information available under separate cover.

#### **Warranty Exclusions:**

DR Drop Protection excluded if applicable. Warranty does not cover damage caused by any act of God, war, terrorism, normal wear and tear, dropping or impact, neglect, power spikes or surges that might potentially damage the equipment's internal components, loss of power, water, smoke, fire, fire suppressant, high winds, lightning strikes, viruses, tampering, operator abuse or use other than what the equipment is intended to be used for. Warranty is non-transferrable. RadSource Imaging Technologies, Inc. is not responsible for lost revenue due to down time.

DR panel drop protection not included in system warranty.

#### **Notes:**

Upon acceptance by both parties, as indicated by their signatures below, this Quotation shall become a binding Sales Contract between BUYER and SELLER.

Installation will be performed during RadSource Imaging Technologies normal operational hours, Monday through Friday, 8:00AM to 5:00PM. Requests for overtime installation srevices will need to be approved by RadSource Imaging Technologies Service Manager. RadSource Imaging Technologies prevailing overtime rates will apply, times the number of men required to perform the required installation.

**Buyer Site Preparation and Installation:** 

Buyer is responsible for all site and/or room preparation including but not limited to Electrical, Plumbing, Overhead Support Work, Power Access, and providing the necessary Networking Access. Professional room drawings can be provided by RadSource Imaging Technologies Inc. upon customer's request with purchase of equipment from RadSource. Drawings will be started once down payment is received. The area where the X-ray equipment is to be installed will need to be completely finished including paint, baseboards, ceiling etc. before the equipment installation can take place.

Quoted Prices and specifications are subject to change without notice prior to signed buyers agreement. Service Agreement Options are available to protect your purchase after the initial warranty has expired.

During or following installation, a RadSource employee may take photographs of the system being installed at your facility along with an image of the outside of your facility to include on their website's installation gallery

(www.radsource.net/installation-gallery/). By signing this agreement, you agree to grant to RadSource, permission to use, reuse, publish, and republish any photographs that they have taken, or may take of your facility, or the equipment installed at your facility. (RadSource will never take photographs of people, or any other hospital, patient or other confidential information.) The photographs can be used in any and all media now or hereafter, including the internet, for the purpose of promoting our installation process at your facility; and to use your company name in connection with any images RadSource so chooses to use. You agree to release and discharge RadSource from any and all claims and demands that may arise out of or in connection with the use of the photographs.

CUSTOMER ACCEPTANCE	RADSOURCE ACCEPTANCE
NAME	NAME
SIGN	SIGN
DATE	DATE

Thank You For Your Business!



8121 NW Terrace, Kansas City, MO 64153 Phone: 816-741-5558 Fax: 816-587-2423 Visit Us Online At: www.radsource.net

### **SALES QUOTATION**

QUOTE # RADQ4231

DATE 11/12/2018

VALID THROUGH 12/12/2018



Account Manager Nancy Briney
Phone 816-741-5558
Email nbriney@radsource.net

### **CUSTOMER DETAILS**

### **Jackson Co. Medical Examiner**

Bill To:

JACKSON CO MED EXAM-2062-MO 950 EAST 21ST STREET Kansas City, MO 64108

Dr Diane Peterson

Phone (816)881-6600

Email KBrooke@jacksongov.org

GPO:

Ship To:

Part Number

Unit Price Qty Ext. Price

New Fuji Prima T2 CR Package Proposal with Cassettes and Imaging Plates

800042146

FCR Prima T2 with FDX Console - Package for PPG Markets

Fujifilm FCR Prima T2 reader unit packaged with the FDX Console workstation and a choice of cassette packages based on specialty market.

.19,000



#### FCR Prima T2 Reader Unit:

- Small Footprint: 3.2 sg. ft, 21"W x 22"D x 15.5"H. 86 lbs.
- Single cassette insertion
- Throughput speeds of up to 47 plates per hour
- 6-outlet surge protection power-strip

#### **Workstation Cart**

- Heavy duty space-saving workstation and countertop workspace, custom designed to fit over the reader making the system a compact all-in-one set up.
- Stitching accessories are available.

### The FDX Console:

Simplifies and speeds workflow with an extra-large image display area, an easy-to-read, customizable interface and an intuitive arrangement of operation buttons to make exams faster for both the technologist and the patient.

- -Hardware includes desktop CPU (Windows 10 Professional, 64-bit, Core i5 processor with 3.2GHz, 4G RAM, and 500 GB HDD), keyboard, mouse, and 19" color touchscreen monitor and barcode scanner.
- -Integrated MC functionality
- -FDX V11.0 Software
- -DICOM CR Store for connectivity to PACS
- -FDX Console Human markers and menus
- -Dynamic Visualization™ advanced image processing presents images with optimal clarity. Includes FNC (Flexible Noise Control) to suppress noise without loss of diagnostic information or sharpness and MFP (Multi-Frequency Processing) providing overall density uniformity for all anatomical regions.
- -Technologist editing tools such as automatic and manual Shutters (black borders) and movable annotation markers
- -Auto trimming simplifies off-center imaging of small anatomy by recognizing the collimated area and applying it to the full screen for optimized display at PACS.
- -QC adjustments including exam reprocessing, sensitivity. latitude, density and contrast
- -Free Text Annotation with commenting text annotation marker capability
- -Image Magnification for enhanced full screen, magnification and zoom image display tools
- -Basic security features include user restriction and are customizable by technologist log in/out.
- -Statistical Analysis reporting of the patient image database, including reason for image rejection coding, for reject or other performance analyses

_		
		•
	-	- 1
٠.	-	_

		System	Purchase Price	\$19,000.00
	Optional Items (Ple	ease add to system pu	ırchase price)	
800019001	14X17 CASSETTE TYPE CC (Option	nal)	\$743.75 × 2	D 148 Initials
800019003	10X12 CASSETTE TYPE CC (Option	nal)	\$481.25 1	Initials
15189368	14 X 17 ST VI IMAGING PLATE (Op	otional)	\$656.25 7 21	De Intias
15189382	10 X 12 ST VI IMAGING PLATE (Op	otional)	\$371.88	Initials
800020606	Prima T Cart (Optional)		\$1,362.50	Initials
Taxes and Freight charges will be included in the final invoice if not			Subtotal	\$19,000.00
specified on t	this quotation.		Sales Tax	
	YES / NO (please circle one		Shipping	
If yes, we wil	Il require a copy of your tax exe	empt certificate.	Grand Total	\$19,000.00
Please choose	e / check by one of the following	g:	Deposit Required	\$0.00
Funding Sour	ce: Cash Credit Card Bank Loan Lease	Credit Card Name: Credit Card Number: Name of Bank: Lessor:		21800
		*If order is being leased, a Le	easing Company's Purchase (	Order is required.

### **Payment Terms:**

Part Number

The terms stated in this agreement are an important part of this business/purchase agreement. The buyer certifies and agrees to the stated payment terms and agrees that they or their purchasing organization will carry them out as stated in this agreement.

50% due with signed buyer's agreement, 40% due 10-business days prior to the scheduled installation of this equipment, remaining 10% due upon completion of installation.

Pricing contained within this proposal excludes shipping and applicable taxes. Shipping charges are pre-paid by RadSource and shall be added to the buyer's final invoice. Any applicable taxes are the responsibility of the buyer. If the buyer is tax-exempt, RadSource Imaging will require a copy of their tax exempt certificate.

System(s) installation will be scheduled on a mutually agreeable date.

#### Warranty:

### RADSOURCE WARRANTY:

Twelve (12) months service from date of completion of installation.

**FUJI WARRANTY:** 

Reader Unit 2 Years
Cassettes/Imaging plates 2 Years
Console 1 Year

All warranties are subject to manufacturer's terms and conditions. Service hours of warranty coverage are Monday through Friday from 8:00 AM to 5:00 PM CST, excluding holidays.

Initial Here

Extended warranty Information available under separate cover.

#### **Warranty Exclusions:**

DR Drop Protection excluded if applicable. Warranty does not cover damage caused by any act of God, war, terrorism, normal wear and tear, dropping or impact, neglect, power spikes or surges that might potentially damage the equipment's internal components, loss of power, water, smoke, fire, fire suppressant, high winds, lightning strikes, viruses, tampering, operator abuse or use other than what the equipment is intended to be used for. Warranty is non-transferrable. RadSource Imaging Technologies, Inc. is not responsible for lost revenue due to down time.

#### Notes:

Upon acceptance by both parties, as indicated by their signatures below, this Quotation shall become a binding Sales Contract between BUYER and SELLER.

Installation will be performed during RadSource Imaging Technologies normal operational hours, Monday through Friday, 8:00AM to 5:00PM. Requests for overtime installation srevices will need to be approved by RadSource Imaging Technologies Service Manager. RadSource Imaging Technologies prevailing overtime rates will apply, times the number of men required to perform the required installation.

#### **Buyer Site Preparation and Installation:**

Buyer is responsible for all site and/or room preparation including but not limited to Electrical, Plumbing, Overhead Support Work, Power Access, and providing the necessary Networking Access. Professional room drawings can be provided by RadSource Imaging Technologies Inc. upon customer's request with purchase of equipment from RadSource. Drawings will be started once down payment is received. The area where the X-ray equipment is to be installed will need to be completely finished including paint, baseboards, ceiling etc. before the equipment installation can take place.

Quoted Prices and specifications are subject to change without notice prior to signed buyers agreement. Service Agreement Options are available to protect your purchase after the initial warranty has expired.

During or following installation, a RadSource employee may take photographs of the system being installed at your facility along with an image of the outside of your facility to include on their website's installation gallery

(www.radsource.net/installation-gallery/). By signing this agreement, you agree to grant to RadSource, permission to use, reuse, publish, and republish any photographs that they have taken, or may take of your facility, or the equipment installed at your facility. (RadSource will never take photographs of people, or any other hospital, patient or other confidential information.) The photographs can be used in any and all media now or hereafter, including the internet, for the purpose of promoting our installation process at your facility; and to use your company name in connection with any images RadSource so chooses to use. You agree to release and discharge RadSource from any and all claims and demands that may arise out of or in connection with the use of the photographs.

CUSTOMER ACCEPTANCE	RADSOURCE ACCEPTANCE
NAME	NAME
SIGN	SIGN
DATE	DATE

Thank You For Your Business!



8121 NW Terrace, Kansas City, MO 64153 Phone: 816-741-5558 Fax: 816-587-2423 Visit Us Online At: www.radsource.net

### **SALES QUOTATION**

QUOTE # RADQ4226

DATE 11/8/2018

VALID THROUGH 12/8/2018



Account Manager Nancy Briney
Phone 816-741-5558
Email nbriney@radsource.net

### **CUSTOMER DETAILS**

### **Jackson County Medical Examine**

Bill To:

JACKSON CO MED EXAM-2062-MO 950 EAST 21ST STREET KANSAS CITY, MO 64108

Dr Diane Peterson / Medical Examin

Ship To:

Phone (816)881-6600

Email KBrooke@jacksongov.org

GPO:

Part Number		Unit Price	Qty	Ext. Price
	New Vieworks VIVIX-S DR 14x17V Wireless	Gadox DR Sy	stem	
HV6056	VIVIX .S View DR 14x17 Wireless Gadox Panel . Scintillator: Amorphous Silicon Gadolinium Oxysulfide Gd.O.S:Tb (Gadox)	\$27,200.00	1	\$27,200.00

Scintillator Amorphous Silicon: Gadolinium Oxysulfide Gd.O.S:Tb (Gadox)

ALTO DR™ 14x17 Wireless Detectors Include: 14x17 Wireless Flat Panel Detector, System Control Unit (SCU), Battery Charger, Battery Pack, Image Processing Viewer Software. Cables, and Acquisition Workstation.

Detection Area: 14x17 inch (35cm x 43cm)

Dimensions: 18.1"W x 15.1"L x .59"D (46 x 38.4 x 1.5 cm)

- Pixels: 2560 x 3072
- Pixel Pitch: 140 im x 140 im
- Scintillator: Cesium (CsI:TI) or Gadox (Gd2O2S:Tb)
- Grayscale: 16,384
- · A/D Conversion: 14 bit
- X-Ray Voltage Range: 40-150 kVp
- Weight: 7.27 lbs./3.3 kg. (CsI) or 6.83 lbs./3.1 kg. (Gadox)
- Power: 7.4V 4000mAh Lithium Ion Polymer Battery



Initial Here

DR ACQUISION WORKSTATION

### **Acquisition Workstation**

Acquisition Hardware:

- Intel Core i3 Processor 6100 CPU (3.7GHz)
- 4 GB, Windows 7 Professional 64-bit
- 1 TB Local Mini-PACs in a RAID 1 (2 drives total)
- Professional Monitor: 22 in. Touch Screen Resolution: 1920 X 1080, full HD Ratio: 1000:1 (typical)

VxVue Imaging Acquisition Software:

- Viewers 5 (LAN) Local Area Network Clinical Viewing Licenses
- AED Technology Automatic Exposure Detection: No Generator Interfacing is required, as the AED Feature instantly detects x-ray upon exposure and then completes the acquisition process for review and transfer of studies to PACS
- DICOM Store/Transfer to PACS, Auto-Routing and Scheduled periodic back-up
- Image Preview in 3 seconds
- Customizable Anatomical exam tree and image processing parameters.
- Image Acquisition + QC (brightness/contrast, window/level, flip, rotate, zoom, annotations, measurements, pan. mag, crop, invert, line measurement, angle measurement, text, arrow, L&R Markers, clockwise rotation, counter-clockwise rotation.)
- Patient Registration Manually or DICOM Modality Worklist. DICOM Modality Worklist is a standard feature and is available for use provided that the EMR/EHR/RIS accommodates Modality Worklist.
- Print: DICOM Printer or Desk Top Paper Printer.
- Stitching: Manual and Auto.
- Archive: CD/DVD Archive to removable drive.
- Export Images as, JPEG, TIFF. BMP, and DICOM with or without Viewer to CD or other media.
- Patient Size Selection: Pediatric, Small, Medium, Large Adult.
- Black Surround Fill

APPS 2.0 2 Days Applications Training Included in purchase price.

**System Purchase Price** 

\$27,200.00

		Unit Price	Qty	Ext. Price
Optional Items (Ple	ease add to system p	urchase price	)	
14X17 WEIGHT DISTRIBUTION CA	P (Optional)	\$400.00	1	Initials
		\$743.75	1	Initials
Please add both prices for the Appropriate Grid set. Specified protective casing.	e complete Protect-A-Grid Prote d Grid is built into a durable pla	ector and	1 _	Initials
TABLETOP DR PANEL HOLDER (Op	otional)	\$600.00	1 _	Initials
(used for performing lateral ta Protect-A-Grid and DR Panel)	bletop radiographs. Holds the			-
ght charges will be included in	the final invoice if not	Subtotal		\$27,200.00
is quotation.		Sales Tax		
		Shipping		
require a copy of your tax exe	empt certificate.	<b>Grand Total</b>		\$27,200.00
/ check by one of the following	g:	Deposit Req	uired	\$0.00
Cash Credit Card Bank Loan Lease	Credit Card Name: Credit Card Number: Name of Bank: Lessor:		28	943.75
	14X17 WEIGHT DISTRIBUTION CA  14X17 PROTECT-A-GRID (PLASTIC HANDLES AND DOUBLE RELEASE  Must be ordered with the below Please add both prices for the Appropriate Grid set. Specified protective casing.  14X17 10:1/215/40-72 GRID (SH TABLETOP DR PANEL HOLDER (O) (used for performing lateral tank Protect-A-Grid and DR Panel)  19th charges will be included in its quotation.  YES / NO (please circle one require a copy of your tax executive casing.  Cash Credit Card Bank Loan	14X17 WEIGHT DISTRIBUTION CAP (Optional)  14X17 PROTECT-A-GRID (PLASTIC) ENCASEMENT WITH 2  HANDLES AND DOUBLE RELEASE LOCKS (Optional)  Must be ordered with the below "Installed Grid" Cat. #K1410  Please add both prices for the complete Protect-A-Grid Protective and properties of the complete Protect-A-Grid Protective casing.  14X17 10:1/215/40-72 GRID (SHORT DIMENSION) (Optional)  TABLETOP DR PANEL HOLDER (Optional)  (used for performing lateral tabletop radiographs. Holds the Protect-A-Grid and DR Panel)  The standard of the final invoice if not is quotation.  YES / NO (please circle one)  require a copy of your tax exempt certificate.  / check by one of the following:	Optional Items (Please add to system purchase price  14X17 WEIGHT DISTRIBUTION CAP (Optional) \$400.00  14X17 PROTECT-A-GRID (PLASTIC) ENCASEMENT WITH 2 \$743.75  HANDLES AND DOUBLE RELEASE LOCKS (Optional)  Must be ordered with the below "Installed Grid" Cat. #K14102040SD.  Please add both prices for the complete Protect-A-Grid Protector and Appropriate Grid set. Specified Grid is built into a durable plastic protective casing.  14X17 10:1/215/40-72 GRID (SHORT DIMENSION) (Optional)  TABLETOP DR PANEL HOLDER (Optional) \$600.00  (used for performing lateral tabletop radiographs. Holds the Protect-A-Grid and DR Panel)  Tight charges will be included in the final invoice if not is quotation.  Subtotal Sales Tax  YES / NO (please circle one) require a copy of your tax exempt certificate.  Check by one of the following:  Check by one of the following:  Credit Card Name:  Credit Card Number:  Bank Loan Name of Bank:	Optional Items (Please add to system purchase price)  14X17 WEIGHT DISTRIBUTION CAP (Optional) \$400.00 1  14X17 PROTECT-A-GRID (PLASTIC) ENCASEMENT WITH 2 \$743.75 1  HANDLES AND DOUBLE RELEASE LOCKS (Optional)  Must be ordered with the below "Installed Grid" Cat. #K14102040SD.  Please add both prices for the complete Protect-A-Grid Protector and Appropriate Grid set. Specified Grid is built into a durable plastic protective casing.  14X17 10:1/215/40-72 GRID (SHORT DIMENSION) (Optional) \$600.00 1  TABLETOP DR PANEL HOLDER (Optional) \$600.00 1  (used for performing lateral tabletop radiographs. Holds the Protect-A-Grid and DR Panel)  Sight charges will be included in the final invoice if not is quotation.  YES / NO (please circle one) require a copy of your tax exempt certificate.  Credit Card Name:  Credit Card Name:  Credit Card Number:  Bank Loan Name of Bank;

### **Payment Terms:**

The terms stated in this agreement are an important part of this business/purchase agreement. The buyer certifies and agrees to the stated payment terms and agrees that they or their purchasing organization will carry them out as stated in this agreement.

50% due with signed buyer's agreement, 40% due 10-Business Days prior to the scheduled installation of this equipment, remaining 10% due upon completion of installation.

Additional convenience fee of 3% will be added to all invoices paid by approved credit card.

Pricing contained within this proposal excludes shipping and applicable taxes. Shipping charges are pre-paid by RadSource and shall be added to the buyer's final invoice. Any applicable taxes are the responsibility of the buyer. If the buyer is tax-exempt, RadSource Imaging will require a copy of their tax exempt certificate.

System(s) installation will be scheduled on a mutually agreeable date. Installation will not begin until the (first and second) progress payments have been submitted to RadSource Imaging Technologies, Inc.

#### Warranty:

#### VIEWORKS ALTO MANUFACTURERS WARRANTY:

Sixty (60) Months on manufacturer defects on the DR Detector. 36-Months on manufacturer defects on the standard PC acquisition station, and software.

#### RADSOURCE LABOR WARRANTY:

Twelve (12) months service from date of completion of installation.

All warranties are subject to manufacturer's terms and conditions. Service hours of warranty coverage are Monday through Friday

RADQ4226

Initial Here

from 8:00 AM to 5:00 PM excluding holidays.

Extended warranty Information available under separate cover.

#### **Warranty Exclusions:**

DR Drop Protection excluded if applicable. Warranty does not cover damage caused by any act of God, war, terrorism, normal wear and tear, dropping or impact, neglect, power spikes or surges that might potentially damage the equipment's internal components, loss of power, water, smoke, fire, fire suppressant, high winds, lightning strikes, viruses, tampering, operator abuse or use other than what the equipment is intended to be used for. Warranty is non-transferrable. RadSource Imaging Technologies, Inc. is not responsible for lost revenue due to down time.

DR panel drop protection not included in system warranty.

#### Notes:

Upon acceptance by both parties, as indicated by their signatures below, this Quotation shall become a binding Sales Contract between BUYER and SELLER.

Installation will be performed during RadSource Imaging Technologies normal operational hours, Monday through Friday, 8:00AM to 5:00PM. Requests for overtime installation srevices will need to be approved by RadSource Imaging Technologies Service Manager. RadSource Imaging Technologies prevailing overtime rates will apply, times the number of men required to perform the required installation.

**Buyer Site Preparation and Installation:** 

Buyer is responsible for all site and/or room preparation including but not limited to Electrical, Plumbing, Overhead Support Work, Power Access, and providing the necessary Networking Access. Professional room drawings can be provided by RadSource Imaging Technologies Inc. upon customer's request with purchase of equipment from RadSource. Drawings will be started once down payment is received. The area where the X-ray equipment is to be installed will need to be completely finished including paint, baseboards, ceiling etc. before the equipment installation can take place.

Quoted Prices and specifications are subject to change without notice prior to signed buyers agreement. Service Agreement Options are available to protect your purchase after the initial warranty has expired.

During or following installation, a RadSource employee may take photographs of the system being installed at your facility along with an image of the outside of your facility to include on their website's installation gallery

(www.radsource.net/installation-gallery/). By signing this agreement, you agree to grant to RadSource, permission to use, reuse, publish, and republish any photographs that they have taken, or may take of your facility, or the equipment installed at your facility. (RadSource will never take photographs of people, or any other hospital, patient or other confidential information.) The photographs can be used in any and all media now or hereafter, including the internet, for the purpose of promoting our installation process at your facility; and to use your company name in connection with any images RadSource so chooses to use. You agree to release and discharge RadSource from any and all claims and demands that may arise out of or in connection with the use of the photographs.

CUSTOMER ACCEPTANCE	RADSOURCE ACCEPTANCE
NAME	NAME
SIGN	SIGN
DATE	DATE

# Thank You For Your Business!